
MASK-OP®

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MASK-OP®, liquid splash resistant and anti-fogging surgical mask, with protective visor, is made of one layer of polyethylene net and three layers of non-woven. These four layers are structured according to the specific function which they are made for: liquid splash resistant, protective, filtering, non-macerating.

MASK-OP® guarantees an effective barrier against the risk of contamination by unexpected splashes of potentially contaminated liquids, which are a possible infection carrier, and so it provides to the wearer the maximum protection against any potential source of contamination. In fact, MASK-OP®'s high resistance to fluids splashes (according to the European norm mentioned below) is guaranteed by a special film, made of a hexagonal structure polyethylene net, which is the second layer of the mask.

The special filter guarantees a "barrier" effect against the biological contamination risk, and also let the wearer breath easily.

ST 1680 – Emission Date – 17th July 2020	Technical development projects: Eric Pirotta		Technical director: Andrea Plebani		Page 1 of 7
Rev. 05 – Rev. Date – 20 Mag. 2021					

The visor for the eyes protection, simple or wrapping, is made of a special polyester film, very transparent, and studied just to avoid any optical distortion. Furthermore, it is submitted to a special treatment that makes the visor anti-fogging on both its internal and external side.

A stripe made of soft expanded polyethylene is bonded at the internal layer of MASK-OP®, on its upper part, in order to provide a soft contact with the face, and also with the function to be like a pillow between the wearer's face and the mask, in order to prevent the warm air inside the mask from going up, obliging the air to go out by the sides of the mask, avoiding the annoying fogging of the glasses and of the visor.

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MASK-OP® is a Type IIR surgical mask, according to the requirements of the norm UNI EN 14683:2019. Furthermore, MASK-OP® is a Class I medical device, self-certificated at the Italian Ministry for Health. It complies with the requirements of UNI EN ISO 10993 (parte 1, 5, 10).

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The Mediberg Disposable Surgical Mask prevents environmental contamination by bacteria and secretions produced by the user and protects it from splashes of liquids. It is intended for use:

- ✓ *by all operating room personnel during surgery, or in other medical activities with similar requirements*
- ✓ *by anyone who needs to reduce the risk of infection spread, in particular in epidemic or pandemic situations*
- ✓ *in the industrial sector, by personnel accessing the cleaning/controlled room or other environments where environmental pollution must be contained.*

ST 1680 – Emission Date – 17th July 2020	Technical development projects: Eric Pirotta		Technical director: Andrea Plebani		Page 2 of 7
Rev. 05 – Rev. Date – 20 Mag. 2021					

COME INDOSSARE LA MASCHERA CHIRURGICA / HOW TO WEAR THE SURGICAL MASK

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It is recommended that the end user inspect the integrity of the product before use. Place the surgical mask on the face making sure to cover the nose and mouth. The opaque, hypoallergenic white side must be facing inwards in contact with the face with the moldable bar at the top. The colored part, also recognizable by the presence of a polished layer, is the outer side of the mask. Place the rubber bands behind the ears and then spread the mask, then tighten the moldable bar around the nose to minimize air spills.

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ST 1680 – Emission Date – 17th July 2020	Technical development projects: Eric Pirotta		Technical director: Andrea Plebani		Page 3 of 7
Rev. 05 – Rev. Date – 20 Mag. 2021					

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ST 1680 – Emission Date – 17th July 2020	Technical development projects: Eric Pirotta		Technical director: Andrea Plebani		Page 4 of 7
Rev. 05 – Rev. Date – 20 Mag. 2021					

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ST 1680 – Emission Date – 17th July 2020	Technical development projects: Eric Pirota		Technical director: Andrea Plebani		Page 5 of 7
Rev. 05 – Rev. Date – 20 Mag. 2021					

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ST 1680 – Emission Date – 17th July 2020	Technical development projects: Eric Pirotta 	Technical director: Andrea Plebani 	Page 6 of 7
Rev. 05 – Rev. Date – 20 Mag. 2021			

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*The mask is distributed by Mediberg S.r.l. and produced in its own production workshop at the following
address: Mediberg S.r.l., Via Vezze, 16/18 - 24050 Calcinate (BG) - Italy*

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*A customer care service is active and it can be reached for whatever information as follows:
Phone: +39 035 83.54.918 – Fax: +39 035 44.29.074 – E-mail: info@mediberg.com*

ST 1680 – Emission Date – 17th July 2020	Technical development projects: Eric Pirotta		Technical director: Andrea Plebani		Page 7 of 7
Rev. 05 – Rev. Date – 20 Mag. 2021					

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**La maschera è distribuita da Mediberg S.r.l. e prodotta nel proprio stabilimento al seguente indirizzo:
Mediberg S.r.l. Via Vezze, 16/18 – 24050 Calcinate (BG) – Italia**